

CLINICAL TOOLS ... CONCLUSIVE RESULTS

The "**R.Test Evolution 4**" from Novacor, is a unique device that has a proven pedigree. It is well established and has been respected throughout the UK since the original R.Test Evolution device was launched in 1995. In fact, previous versions of the R.Test Evolution are currently in daily use in more than 300 NHS and Private Hospitals and an ever growing number of GP Surgeries & Primary Care Centers throughout the UK and many more around the world.

The "**R.Test Evolution 4**" takes less than a minute to set-up and is extremely easy to fit to your patient. It is exceptionally light and comfortable to wear and most patients genuinely forget that they are wearing it at all.

The "**R.Test Evolution 4**" is the only truly **AUTOMATIC**, Non-Invasive, ECG Arrhythmia Detection Device. This means that it can be used successfully to screen for infrequent and/or asymptomatic arrhythmia's, including those that occur during sleep. It can be fitted to virtually any patient in the knowledge that results are not dependant upon the patient physically activating the device.... during symptoms or otherwise. Of course , patients can mark symptomatic events if they wish and are able to.

The "**R.Test Evolution 4**" can be worn from the traditional 24-hour period or less, or now for anything up to **32 days!** whilst continuously monitoring, detecting and storing abnormal ECG. Unsurprisingly, this extended monitoring period is scientifically proven to yield far more "conclusive" results than the traditional 24-hour tape test, and an excellent step before using very expensive, single use, invasive implantable loop recorders, with their associated procedures, complications and costs.